

FP 4 MINI



Fisiopress "mini" is a sequential, simple, economic and effective home care compression therapy.

It is within the reach of all those who need it.


THE MINI SYSTEM HAS A HOMOGENEOUSLY MODULATED PRESSURE GRADIENT, IN ORDER TO EXERCISE A CONSTANT PRESSURE ON THE TREATED LIMB, FROM THE PERIPHERY TO THE BASE. REGARDLESS OF ITS GEOMETRY (VOLUME). THIS GUARANTEES A TRUE EFFECTIVE AND PLEASANT DRAINAGE.

The sequential compression therapy at home represents an effective treatment for problems of edema, lymphedema, swollen legs, cellulite, and perform a correct venous and lymphatic circulation by improving muscle tension related to stress or pain: its action of compression and decompression of the pneumatic sectors normalizes and restores the venous and lymphatic return and increases tissue spraying.

USE AND APPLICATION

The Fisiopress FP 4 MINI device is the ideal device for:

- venous insufficiency and its consequences
- venous edema
- forms of stagnation of liquids in the lower and upper limbs
-) Water retention
- particularly indicated also in post traumatic oedemas, ulcers
-) and in the post course operative
- after ligation of the perforators or stripping of the saphenous vein
- lymphatic edema and lymphedema
- cellulite
- aching and heavy legs
- recovery of muscle tone in immobilized or disabled patients



TECHNICAL FEATURES

-) Distal to proximal pressure gradient sequential pressotherapy.
-) Portable-professional address system, equipped with 4 sequential outputs, capable of treat one or two limbs simultaneously.
- Furthermore, it is possible to contemporary connect a pair of leg sleeves and an abdominal strip
-) Pressure range: 0 - 80 mm of Hg.

-) Operating cycle: truly sequential.
-) DASC system (Distribution to communicating sectors): this system, the only Fisiopress, allows you to have a homogeneous pressure in all the pneumatic chambers as you go fill with air, thus allowing the garments to adapt to the geometry of the **treated limb**, compressing it uniformly without danger of stagnation or reverse flows.
-) Entirely hand-sewn garments.
-) Garment available for the lower limb: leg sleeve (4 sectors)
-) Abdominal strip (2 sectors)
-) Abdominal band (4 sectors) to be used individually.
-) Extender to increase the circular dimension of the garment.
-) Garment available for upper limb: arm sleeve (4 sectors).
-) All garments have parallel and partially overlapped internal sectors.
-) Possibility of individual replacement of each internal pneumatic bag.
-) Entirely hand-sewn garments.
-) Medical CE mark





FP MINI WITH DOUBLE LEG SLEEVES



FP MINI WITH DOUBLE LEG SLEEVES AND ABDOMINAL STRIP



FP 4 MINI WITH SINGLE ARM SLEEVE