

fisiopress

FP 4/8

COMPRESSION THERAPY



Professional devices for home care use represent a new segment for pressure therapy equipment as they represent the perfect combination of professional appliance and home appliance. The goal is to bring the characteristics of the best professional systems to the home setting, thus giving the user a truly effective tool for their needs. The application of sequential compression therapy represents an effective treatment for the problems of edema, lymphedema, swollen legs, cellulite, able to promote correct venous and lymphatic circulation by improving muscle tension linked to stress or pain: its compression and decompression of the pneumatic sectors normalizes and restores venous and lymphatic return and increases tissue circulation.

USE AND INDICATIONS

This equipment has been designed to address, at home, important pathologies that require pressure therapy systems of more than 4 or 6 sectors. The FP 4-8 device with its 8 compression sectors is a sequential home system with features of functionality and effectiveness derived from professional systems.

The Fisiopress FP 4/8 system is the ideal device for:

- venous pathologies (venous insufficiency and its consequences)
- venous edema
- lymphatic pathologies
- lymphatic edema
- secondary lymphedema of the upper limb
- in all forms of fluid stagnation of the lower and upper limbs,
- post traumatic edema, in ulcers and in the postoperative course
- after ligation of the perforators or stripping of the saphenous vein
- cellulitis
- painful and heavy legs
- recovery of muscle tone in immobilized or disabled patients

TECHNICAL FEATURES

- Sequential pressure gradient distal-proximal pressure therapy.
- Portable-professional address system, equipped with 8 sequential outputs, capable of treating one or two limbs at the same time.
- Possibility, moreover, to connect a pair of leggings and the abdominal belt at the same time.
- Pressure range: 0 - 160 mm of Hg. on one garment.
- Operating cycle: 60 seconds.
- D.A.S.C. system (Distribution to communicating sectors): this system, the only Fisiopress, allows to have an always homogeneous pressure in all the chambers as they fill with air, thus allowing the terminals to adapt to the geometry of the treated limb and compress it uniformly without danger of stagnation or reverse flows.
- Garments available for lower limb: 8 sectors leg.
- Abdominal band (2 sectors)
- Expander to increase the circular size of the leg.
- Garment available for upper limb: cuff (8 sectors).
- All garments have parallel and partially overlapping internal sectors.
- Possibility of individual replacement of each internal pneumatic bag.
- CE medical mark